

Aprons

■ Characteristics

- + Good fit and high comfort are insured by broad straps
- + Easy cleaning without residue
- + Aprons are detectable
- + Blue hygienic TPU-X-straps
- + If pointed knives are used, type 4.2 is highly recommended!

Wire gage:
7,0 mm or 4,2 mm

DIN EN ISO 13998:2003
(performance level 2)



Type	Measure	Weight	Height	Chest Measure	Chainmail dimensions	Item code
Ø 7.0 x 0.7	45 x 75 cm	1,1 kg	160 – 172 cm	80 – 100 cm	7,0 x 0,7 mm	AS 02 4575 0000
Ø 7.0 x 0.7	55 x 80 cm	1,3 kg	174 – 186 cm	100 – 115 cm	7.0 x 0,7 mm	AS 02 5580 0000
Ø 7.0 x 0.7	60 x 85 cm	1,45 kg	188 – 196 cm	120 – 130 cm	7,0 x 0,7 mm	AS 02 6085 0000
Ø 4.2 x 0.55	45 x 75 cm	1,3 kg	160 – 172 cm	80 – 100 cm	4,2 x 0,55 mm	AS 10 4575 0000
Ø 4.2 x 0.55	55 x 80 cm	1,5 kg	174 – 186 cm	100 – 115 cm	4,2 x 0,55 mm	AS 10 5580 0000
Ø 4.2 x 0.55	60 x 85 cm	1,7 kg	188 – 196 cm	120 – 130 cm	4,2 x 0,55 mm	AS 10 6085 0000

Notes for cleaning and storage:

After use, aprons and straps are to be cleaned with hot water 50° C with a standard detergent commonly used in the meat and food industry. We recommend cleaning with a standard washing machine (laundry). After cleaning, rinse the aprons with 85° C hot water. After drying, the aprons should be sprayed with a disinfectant, allowing the disinfectant to settle on the apron for several minutes. After cleaning aprons might be dried and stored in a dry, cool and well ventiated place.



niroflex aprons also
available in a
detectable version



DIN EN ISO
13998:2003

